

Anti-Stk39 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Uew8

Item Number: CM0006546



Introduction

STK39 Rabbit mAb (Q9UEW8) validated for WB, IF/ICC, IHC-P, ELISA. High specificity for human, mouse, rat STK39 kinase, a WNK-SPAK/OSR1 cascade component regulating ion transport and blood pressure.

[Learn More](#)

Parameter	Value
Product Name	STK39 Rabbit mAb
Remarks/Alias	DCHT; PASK; SPAK; Ste 20 related kinase; STK39
Species	Human
GeneID (Human)	27347
GeneID	27347
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 401-545 of human STK39 (Q9UEW8).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UEW8
Protein Weight	59kDa
Shipping	Ice bag