

Anti-Stk39 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uew8

Item Number: CM0026944



Introduction

Rabbit polyclonal antibody against human STK39 kinase, validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for studying the WNK-SPAK/OSR1 kinase cascade in ion transport and blood pressure regulation.

[Learn More](#)

Parameter	Value
Product Name	STK39 Rabbit pAb
Remarks/Alias	DCHT; PASK; SPAK; STK39
Species	Human
GeneID (Human)	27347
GeneID	27347
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human STK39 (Q9UEW8).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UEW8
Protein Weight	59kDa
Shipping	Ice bag