

Anti-Stk3/Mst2 Monoclonal Antibody For Wb, If/Icc, Elisa - Q13188

Item Number: CM0010094



Introduction

High-quality rabbit monoclonal antibody against human STK3/MST2, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for apoptosis and Hippo signaling research.

[Learn More](#)

Parameter	Value
Product Name	STK3/MST2 Rabbit mAb
Remarks/Alias	KRS1, MST2, STK3/MST2
Species	Human
Gene ID (Human)	6788
Gene ID	6788
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human STK3/MST2 (Q13188).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13188
Protein Weight	56kDa
Shipping	Ice bag