

Anti-Stk3/Mst2 (Phospho T117) Polyclonal Antibody For Wb And Elisa - Q13188

Item Number: CM0001399



Introduction

Rabbit polyclonal antibody targeting phosphorylated STK3/MST2 at Thr117, validated for Western blot and ELISA. Detects human, mouse, and rat STK3/MST2. Ideal for apoptosis and Hippo signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STK3/MST2-T117 Rabbit pAb
Remarks/Alias	KRS1; MST2; Phospho-STK3/MST2-T117
Species	Human
Gene ID (Human)	6788
Gene ID	6788
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13188
Protein Weight	56kDa
Shipping	Ice bag