

Anti-Stk3/Mst2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13188

Item Number: CM0028448



Introduction

Rabbit polyclonal antibody targeting human STK3/MST2, suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying Hippo pathway, apoptosis, and centrosome function.

[Learn More](#)

Parameter	Value
Product Name	STK3/MST2 Rabbit pAb
Remarks/Alias	KRS1; MST2; STK3/MST2
Species	Human
GeneID (Human)	6788
GeneID	6788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 323-450 of human STK3/MST2 (NP_006272.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13188
Protein Weight	56kDa
Shipping	Ice bag