

Anti-Stmn2 Rabbit Polyclonal Antibody For Wb, Elisa - Q93045

Item Number: CM0030443

Introduction

Rabbit polyclonal antibody against STMN2 (Stathmin-2/SCG10). Validated in WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for neuroscience research on microtubule dynamics and neuronal development.

[Learn More](#)

Parameter	Specification
Product Name	STMN2 Rabbit pAb
Remarks/Alias	SCG10; SCGN10; STMN2
Species	Human
GeneID (Human)	11075
GeneID	11075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human STMN2 (NP_008960.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93045
Protein Weight	21kDa
Shipping	Ice bag