


Anti-Stra13 Polyclonal Antibody For Wb, Elisa - A8Mt69

Item Number: CM0030306



Introduction

Rabbit anti-human STRA13/CENPX polyclonal antibody validated for WB and ELISA. Suitable for research on Fanconi anemia, kinetochore assembly, and DNA repair. — A8MT69

[Learn More](#)

Parameter	Value
Product Name	STRA13 Rabbit pAb
Remarks/Alias	D9, MHF2, CENP-X, FAAP10, STRA13
Species	Human
Gene ID	201254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human STRA13 (NP_659435.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	A8MT69
Protein Weight	9kDa
Shipping	Ice bag