


## Anti-Strada Rabbit Polyclonal Antibody For Wb And Elisa - Q7Rtn6

Item Number: CM0023095



### Introduction

Rabbit polyclonal antibody against human STRADA (STE20-related kinase adapter protein alpha). Validated for WB and ELISA, with mouse cross-reactivity. Suitable for studying LKB1/STK11 activation and cell cycle regulation.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | STRADA Rabbit pAb                                                                                                  |
| Remarks/Alias    | LYK5; PMSE; Stik; STRAD; NY-BR-96; STRAD alpha; STRADA                                                             |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 92335                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 25-210 of human STRADA (NP_699166.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Mouse                                                                                                              |
| SWISS            | Q7RTN6                                                                                                             |
| Protein Weight   | 48kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |