

## Anti-Strap Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y3F4

Item Number: CM0025237

### Introduction

High-quality rabbit polyclonal antibody against human STRAP protein. Validated for WB, IF/ICC, ELISA applications. Detects endogenous STRAP in human, mouse, and rat samples. Suitable for spliceosome and TGF-beta signaling pathway research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | STRAP Rabbit pAb                                                                                                 |
| Remarks/Alias    | MAWD; PT-WD; UNRIP; STRAP                                                                                        |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 11171                                                                                                            |
| GeneID           | 11171                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human STRAP (NP_009109.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q9Y3F4                                                                                                           |
| Protein Weight   | 38kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |