

Anti-Strap Rabbit Polyclonal Antibody For Wb - Q9Y3F4

Item Number: CM0023283

Introduction

Rabbit polyclonal antibody against human STRAP (MAWD/UNRIP), validated by knockout, for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Detects a 38 kDa band. Ideal for studies on snRNP assembly and TGF-beta signaling.

[Learn More](#)

Parameter	Description
Product Name	[KO Validated] STRAP Rabbit pAb
Remarks/Alias	MAWD; PT-WD; UNRIP; AP
Species	Human
GeneID (Human)	11171
GeneID	11171
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human STRAP (NP_009109.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3F4
Protein Weight	38kDa
Shipping	Ice bag