


Anti-Strbp Polyclonal Antibody For Wb And Elisa - Q96Si9

Item Number: CM0014040



Introduction

Anti-STRBP rabbit polyclonal antibody validated for WB and ELISA. Targets human spermatid perinuclear RNA-binding protein, involved in spermatogenesis and RNA binding.

[Learn More](#)

Parameter	Value
Product Name	STRBP Rabbit pAb
Remarks/Alias	p74; SPNR; ILF3L; HEL162; STRBP
Species	Human
GeneID (human)	55342
GeneID	55342
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-190 of human STRBP (NP_060857.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96SI9
Protein Weight	74kDa
Shipping	Ice bag