


Anti-Strip1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Vsl9

Item Number: CM0018179



Introduction

Anti-STRIP1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA with human, mouse, rat reactivity. Research use in cytoskeletal organization and STRIPAK complex studies targeting UniProt-reviewed Q5VSL9 protein.

[Learn More](#)

Parameter	Value
Product Name	STRIP1 Rabbit pAb
Remarks/Alias	FAM40A; FAR11A; STRIP1
Species	Human
GeneID (Human)	85369
GeneID	85369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 743-837 of human STRIP1 (NP_149079.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VSL9
Protein Weight	96kDa
Shipping	Ice bag