

Anti-Strn Rabbit Polyclonal Antibody For Wb - 043815

Item Number: CM0029025

Introduction

Rabbit polyclonal antibody against human STRN, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying STRIPAK complex signaling in cell growth and differentiation.

[Learn More](#)

Parameter	Specification
Product Name	STRN Rabbit pAb
Remarks/Alias	SG2NA; STRN1; PPP2R6A; STRN
Species	Human
GeneID (human)	6801
GeneID	6801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 581-780 of human STRN (NP_003153.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43815
Protein Weight	86kDa
Shipping	Ice bag