


Anti-Strn3 Polyclonal Antibody For Wb - Q13033

Item Number: CM0028802



Introduction

High-quality rabbit polyclonal antibody against human STRN3 (SG2NA), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for STRIPAK complex and signaling research.

[Learn More](#)

Parameter	Value
Product Name	STRN3 Rabbit pAb
Remarks/Alias	SG2NA; S/G2NA; PPP2R6B; STRN3
Species	Human
GeneID	29966
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-460 of human STRN3 (NP_001077362.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13033
Protein Weight	87kDa
Shipping	Ice bag