

Anti-Stt3B Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tcj2

Item Number: CM0013725



Introduction

Rabbit polyclonal antibody targeting human STT3B, the catalytic subunit of the OST complex, validated for WB and ELISA. Detects human and mouse STT3B (94 kDa). Ideal for N-glycosylation and ERAD research.

[Learn More](#)

Parameter	Value
Product Name	STT3B Rabbit pAb
Remarks/Alias	SIMP; CDG1X; STT3-B; STT3B
Species	Human
Gene ID (Human)	201595
Gene ID	201595
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human STT3B (NP_849193.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TCJ2
Protein Weight	94kDa
Shipping	Ice bag