


Anti-Stub1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Une7

Item Number: CM0014938



Introduction

CamelBio's STUB1 rabbit polyclonal antibody, validated for WB, IF/ICC, IP, and ELISA in human, mouse, and rat samples. Targets the E3 ubiquitin ligase CHIP involved in protein quality control and mitochondrial homeostasis.

[Learn More](#)

Parameter	Value
Product Name	STUB1 Rabbit pAb
Remarks/Alias	CHIP; SCA48; UBOX1; SCAR16; HSPABP2; NY-CO-7; SDCCAG7; STUB1
Species	Human
GeneID (Human)	10273
GeneID	10273
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human STUB1 (NP_005852.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNE7
Protein Weight	35kDa
Shipping	Ice bag