

Anti-Stx12 Rabbit Polyclonal Antibody For Wb And Elisa - Q86Y82

Item Number: CM0028148



Introduction

Anti-STX12 Rabbit Polyclonal Antibody validated for Western Blot (WB) and ELISA. Reacts with Human, Mouse, Rat. Targets syntaxin-12 involved in endocytic recycling and AMPAR trafficking. Ideal for neuroscience and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	STX12 Rabbit pAb
Remarks/Alias	STX13; STX14; STX12
Species	Human
GeneID (Human)	23673
GeneID	23673
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-276 of human STX12 (NP_803173.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86Y82
Protein Weight	32kDa
Shipping	Ice bag