

Anti-Stx1A Polyclonal Antibody For Wb, If/Icc, Elisa - Q16623

Item Number: CM0029231

Introduction

Rabbit polyclonal antibody against human STX1A (Syntaxin-1A, HPC-1), validated for WB, IF/ICC, ELISA. Cross-reacts with human and rat. Ideal for studying SNARE-mediated exocytosis.

[Learn More](#)

Parameter	Value
Product Name	STX1A Rabbit pAb
Remarks/Alias	STX1; HPC-1; P35-1; SYN1A; STX1A
Species	Human
GeneID (human)	6804
GeneID	6804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-288 of human STX1A (NP_004594.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q16623
Protein Weight	33kDa
Shipping	Ice bag