

## Anti-Stx2 Rabbit Polyclonal Antibody For Wb, Elisa - P32856

Item Number: CM0018317

### Introduction

STX2 Rabbit Polyclonal Antibody validated for WB, ELISA applications with cross-reactivity to Human, Mouse, Rat. Ideal for epithelial morphogenesis and exocytosis research. UniProt P32856.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | STX2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | EPM; EPIM; STX2A; STX2B; STX2C; STX2                                                                            |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 2054                                                                                                            |
| Gene ID          | 2054                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human STX2 (NP_001971.2) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P32856                                                                                                          |
| Protein Weight   | 33kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |