

Anti-Stx6 Rabbit Polyclonal Antibody For Wb, Elisa - 043752

Item Number: CM0022227



Introduction

High-quality rabbit polyclonal antibody targeting human Syntaxin-6 (STX6), validated for Western blot and ELISA applications. Exhibits cross-reactivity with mouse. Suitable for investigating Golgi apparatus, endosome trafficking, and membrane fusion mechanisms.

[Learn More](#)

Parameter	Value
Product Name	STX6 Rabbit pAb
Remarks/Alias	STX6
Species	Human
GeneID (Human)	10228
GeneID	10228
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human STX6 (NP_005810.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43752
Protein Weight	29kDa
Shipping	Ice bag