

Anti-Stx8 Polyclonal Antibody For Wb, Elisa - Q9Unk0

Item Number: CM0018485

Introduction

Anti-STX8 polyclonal antibody validated for Western blot (WB) and ELISA, recognizing human, mouse, and rat STX8. Ideal for studying retrograde transport from cis-Golgi to ER. Rabbit-derived recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	STX8 Rabbit pAb
Remarks/Alias	CARB; STX8
Species	Human
GeneID (Human)	9482
GeneID	9482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human STX8 (NP_004844.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNK0
Protein Weight	27kDa
Shipping	Ice bag