

Anti-Stxbp1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P61764

Item Number: CM0009838



Introduction

High-quality STXBP1 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. 68 kDa protein, ideal for synaptic function and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	STXBP1 Rabbit mAb
Remarks/Alias	P67; DEE4; NSEC1; UNC18; N-Sec1; RBSEC1; unc-18A; unc18-1; MUNC18-1
Species	Human
GeneID (Human)	6812
GeneID	6812
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61764
Protein Weight	68 kDa
Shipping	Ice bag