

Anti-Stxbp1 Polyclonal Antibody For If/Icc, If-P, Elisa - P61764

Item Number: CM0026136



Introduction

High-quality STXBP1 rabbit polyclonal antibody validated for IF/ICC, IF-P, and ELISA. Detects human and rat STXBP1 (MUNC18-1), a key synaptic vesicle fusion regulator. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	STXBP1 Rabbit pAb
Remarks/Alias	P67; DEE4; NSEC1; UNC18; N-Sec1; RBSEC1; unc-18A; unc18-1; MUNC18-1; STXBP1
Species	Human
Gene ID (Human)	6812
Gene ID	6812
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P61764
Protein Weight	68kDa
Shipping	Ice bag