


Anti-Stxbp2 Rabbit Polyclonal Antibody For Wb And Elisa - Q15833

Item Number: CM0029645



Introduction

Rabbit polyclonal antibody targeting STXBP2 (Munc18-2), validated for Western blot (WB) and ELISA in human, mouse, and rat samples. Ideal for research on SNARE-mediated vesicle trafficking, exocytosis, and immune cell cytotoxic granule biology.

[Learn More](#)

Parameter	Value
Product Name	STXBP2 Rabbit pAb
Remarks/Alias	FHL5; UNC18B; Hunc18b; UNC18-2; pp10122; unc-18B; MUNC18-2; STXBP2
Species	Human
GeneID (Human)	6813
GeneID	6813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human STXBP2 (NP_001120868.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15833
Protein Weight	66kDa
Shipping	Ice bag