


Anti-Stxbp2 Rabbit Polyclonal Antibody For Wb And Elisa - Q15833

Item Number: CM0014465



Introduction

Rabbit polyclonal antibody targeting human STXBP2 (Unc-18 homolog 2, Munc18-2), suitable for Western blotting and ELISA. Cross-reacts with human and mouse STXBP2, supporting intracellular vesicle trafficking and immune cell exocytosis studies.

[Learn More](#)

Parameter	Value
Product Name	STXBP2 Rabbit pAb
Remarks/Alias	FHL5; UNC18B; Hunc18b; UNC18-2; pp10122; unc-18B; MUNC18-2; STXBP2
Species	Human
GeneID (Human)	6813
GeneID	6813
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15833
Protein Weight	66kDa
Shipping	Ice bag