


Anti-Stxbp4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Zwj1

Item Number: CM0027243



Introduction

Anti-STXBP4 polyclonal antibody validated for WB, IHC-P, and ELISA in human and mouse. Targets syntaxin-binding protein 4 (Synip), a key regulator of vesicle trafficking and insulin signaling. Ideal for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	STXBP4 Rabbit pAb
Remarks/Alias	Synip; STXBP4
Species	Human
GeneID (Human)	252983
GeneID	252983
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 344-553 of human STXBP4 (NP_848604.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6ZWJ1
Protein Weight	62kDa
Shipping	Ice bag