


Anti-Styx Polyclonal Antibody For Wb, Elisa - Q8Wuj0

Item Number: CM0015175



Introduction

Rabbit polyclonal antibody targeting human STYX (25 kDa), suitable for WB and ELISA. Ideal for investigating MAPK/ERK signaling, FBXW7-mediated degradation, and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	STYX Rabbit pAb
Remarks/Alias	STYX
Species	Human
GeneID (Human)	6815
GeneID	6815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human STYX (NP_660294.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WUJ0
Protein Weight	25kDa
Shipping	Ice bag