

Anti-Sub1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P53999

Item Number: CM0029929



Introduction

Anti-SUB1 rabbit polyclonal antibody detects the transcriptional coactivator p15 (PC4) with applications in WB, IF/ICC, IP, and ELISA. Cross-reacts with human and mouse. Expected protein weight 14 kDa.

[Learn More](#)

Parameter	Value
Product Name	SUB1 Rabbit pAb
Remarks/Alias	P15; PC4; p14; SUB1
Species	Human
Gene ID (Human)	10923
Gene ID	10923
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human SUB1 (NP_006704.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P53999
Protein Weight	14kDa
Shipping	Ice bag