

Anti-Sucla2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P2R7

Item Number: CM0018209



Introduction

High-quality anti-SUCLA2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets mitochondrial succinate-CoA ligase subunit beta, crucial for TCA cycle and ATP synthesis.

[Learn More](#)

Parameter	Value
Product Name	SUCLA2 Rabbit pAb
Remarks/Alias	A-SCS; A-BETA; MTDPS5; LINC00444; SCS-betaA; SUCLA2
Species	Human
GeneID (human)	8803
GeneID	8803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SUCLA2 (NP_003841.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2R7
Protein Weight	50kDa
Shipping	Ice bag