

Anti-Suc1g2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96I99

Item Number: CM0025380



Introduction

KO validated anti-SUCLG2 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat SUCLG2 (GTP-specific succinyl-CoA synthetase), a key TCA cycle enzyme. Suited for mitochondrial metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SUCLG2 Rabbit pAb
Remarks/Alias	G-SCS; GBETA; GTPSCS; G2
Species	Human
Gene ID (Human)	8801
Gene ID	8801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-300 of human SUCLG2 (NP_003839.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96I99
Protein Weight	47kDa
Shipping	Ice bag