

Anti-Sugct Polyclonal Antibody For Wb, Elisa - Q9Hac7

Item Number: CM0026100



Introduction

SUGCT Rabbit pAb is a polyclonal antibody validated for WB and ELISA. It reacts with human, mouse, and rat SUGCT. The target is a mitochondrial CoA-transferase involved in lysine/tryptophan degradation. Ideal for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SUGCT Rabbit pAb
Remarks/Alias	GA3; ORF19; DERP13; C7orf10; SUGCT
Species	Human
GeneID (human)	79783
GeneID	79783
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human SUGCT (NP_001180240.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HAC7
Protein Weight	48kDa
Shipping	Ice bag