


Anti-Sugp2 Polyclonal Antibody For Wb, Elisa - Q8Ix01

Item Number: CM0018507



Introduction

SUGP2 Rabbit pAb is a polyclonal antibody for WB and ELISA applications, recognizing human and mouse SUGP2. This splicing factor localizes to the nucleus and is expressed in testis, brain, and kidney. Ideal for splicing research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SUGP2 Rabbit pAb
Remarks/Alias	SFRS14; SRSF14; SUGP2
Species	Human
GeneID (Human)	10147
GeneID	10147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human SUGP2 (NP_055699.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IX01
Protein Weight	120kDa
Shipping	Ice bag