

Anti-Sulf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Iwu6

Item Number: CM0014049

Introduction

Rabbit polyclonal antibody targeting human SULF1, validated for Western blot and ELISA. Detects SULF1 in human and rat samples. Ideal for apoptosis and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	SULF1 Rabbit pAb
Remarks/Alias	SULF-1; SULF1
Species	Human
GeneID (human)	23213
GeneID	23213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-610 of human SULF1 (NP_001121676.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IWU6
Protein Weight	101kDa
Shipping	Ice bag