

Anti-Sult1A1 Rabbit Polyclonal Antibody For Wb And Elisa - P50225

Item Number: CM0020325



Introduction

Anti-SULT1A1 rabbit polyclonal antibody validated for Western blotting and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying xenobiotic and hormone metabolism.

[Learn More](#)

Parameter	Value
Product Name	SULT1A1 Rabbit pAb
Remarks/Alias	PST; STP; STP1; P-PST; ST1A1; ST1A3; TSPST1; ts-PST; P-PST 1; HAST1/HAST2; SULT1A1
Species	Human
GeneID (human)	6817
GeneID	6817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human SULT1A1 (NP_803565.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50225
Protein Weight	34kDa
Shipping	Ice bag