

Anti-SULT1C4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075897

Item Number: CM0021898



Introduction

Rabbit polyclonal antibody against human SULT1C4. Suitable for Western blot, immunofluorescence, and ELISA. Detects endogenous levels of SULT1C4 in human, mouse, and rat samples. Targets a sulfotransferase involved in xenobiotic and drug metabolism.

[Learn More](#)

Parameter	Value
Product Name	SULT1C4 Rabbit pAb
Remarks/Alias	SULT1C; SULT1C2; SULT1C4
Species	Human
Gene ID	27233 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human SULT1C4 (NP_006579.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075897
Protein Weight	36kDa
Shipping	Ice bag