

Anti-Sult1E1 Polyclonal Antibody For Wb, If/Icc, Elisa - P49888

Item Number: CM0028292



Introduction

High-quality anti-SULT1E1 polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Targets human SULT1E1, with cross-reactivity to mouse and rat. Ideal for estrogen metabolism and hormone-related research.

[Learn More](#)

Parameter	Value
Product Name	SULT1E1 Rabbit pAb
Remarks/Alias	EST; STE; EST-1; ST1E1; SULT1E1
Species	Human
GeneID (Human)	6783
GeneID	6783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-294 of human SULT1E1 (NP_005411.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49888
Protein Weight	35kDa
Shipping	Ice bag