


## Anti-Sult2A1 Polyclonal Antibody For Wb, Elisa - Q06520

Item Number: CM0027745



### Introduction

SULT2A1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse. Detects ~34 kDa endogenous SULT2A1, a key steroid and bile acid sulfotransferase.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SULT2A1 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST; SULT2A3; DHEA-ST8; SULT2A1                                                           |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 6822                                                                                                                                    |
| GeneID           | 6822                                                                                                                                    |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human SULT2A1 (NP_003158.2). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, ELISA                                                                                                                               |
| Cross Reactivity | Human, Mouse                                                                                                                            |
| SWISS            | Q06520                                                                                                                                  |
| Protein Weight   | 34kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |