


Anti-Sult2A1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q06520

Item Number: CM0005187



Introduction

SULT2A1 Rabbit Monoclonal Antibody validated for WB, IF-P, ELISA. Reacts with human and rat SULT2A1 (34kDa), a key steroid and bile acid sulfotransferase. Ideal for metabolism and detoxification studies.

[Learn More](#)

Parameter	Value
Product Name	SULT2A1 Rabbit mAb
Remarks/Alias	HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST; SULT2A3; DHEA-ST8; SULT2A1
Species	Human
Gene ID (Human)	6822
Gene ID	6822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-285 of human SULT2A1 (Q06520)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q06520
Protein Weight	34kDa
Shipping	Ice bag