


Anti-Sult2A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q06520

Item Number: CM0016662



Introduction

Anti-SULT2A1 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat SULT2A1; immunogen recombinant protein (aa1-145). Suitable for steroid sulfonation and bile acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SULT2A1 Rabbit pAb
Remarks/Alias	HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST; SULT2A3; DHEA-ST8; SULT2A1
Species	Human
GeneID (human)	6822
GeneID	6822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human SULT2A1 (NP_003158.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06520
Protein Weight	34kDa
Shipping	Ice bag