

## Anti-Sult2B1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00204

Item Number: CM0028426



### Introduction

CamelBio's Anti-SULT2B1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA, reactive to human, mouse, rat. Targets sulfotransferase 2B1 involved in steroid sulfation and cholesterol metabolism. Ideal for lipid metabolism research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | SULT2B1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | HSST2; ARCI14; SULT2B1                                                                                             |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 6820                                                                                                               |
| Gene ID          | 6820                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-365 of human SULT2B1 (NP_814444.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O00204                                                                                                             |
| Protein Weight   | 41kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |