


Anti-Sumf2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nbj7

Item Number: CM0021188



Introduction

Rabbit polyclonal antibody against SUMF2 (pFGE), an inactive sulfatase-modifying factor that inhibits SUMF1. Suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Immunogen: recombinant protein (aa 45-320).

[Learn More](#)

Parameter	Value
Product Name	SUMF2 Rabbit pAb
Remarks/Alias	pFGE; SUMF2
Species	Human
GeneID (Human)	25870
GeneID	25870
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-320 of human SUMF2 (NP_056226.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBJ7
Protein Weight	34kDa
Shipping	Ice bag