


Anti-Sumo1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P63165

Item Number: CM0019996



Introduction

Anti-SUMO1 Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, ELISA. Detects SUMO1 in Human, Mouse, Rat. Ubiquitin-like modifier involved in nuclear transport, DNA repair, and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	SUMO1 Rabbit pAb
Remarks/Alias	DAP1; GMP1; PIC1; SMT3; UBL1; OFC10; SENP2; SMT3C; SMT3H3; SUMO1
Species	Human
GeneID (Human)	7341
GeneID	7341
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human SUMO1 (NP_003343.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63165
Protein Weight	12kDa
Shipping	Ice bag