

Anti-Sumo1 Rabbit Monoclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Elisa - P63165

Item Number: CM0002520

Introduction

Rabbit monoclonal antibody targeting SUMO1 (P63165), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential for studying SUMOylation in nuclear transport, DNA repair, and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SUMO1 Rabbit mAb
Remarks/Alias	DAP1; GMP1; PIC1; SMT3; UBL1; OFC10; SENP2; SMT3C; SMT3H3; O1
Species	Human
GeneID	7341 (Human)
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-101 of human Sumo 1 (P63165).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63165
Protein Weight	12kDa
Shipping	Ice bag