


Anti-Sumo1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - P63165

Item Number: CM0023547



Introduction

KO-validated polyclonal antibody targeting human SUMO1, suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying SUMOylation in nuclear transport, DNA repair, and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SUMO1 Rabbit pAb
Remarks/Alias	DAP1, GMP1, PIC1, SMT3, UBL1, OFC10, SENP2, SMT3C, SMT3H3, O1
Species	Human
Gene ID	7341
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SUMO1 (NP_003343.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63165
Protein Weight	12kDa
Shipping	Ice bag