

Anti-Sumo2/3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P61956/P55854

Item Number: CM0016585



Introduction

Rabbit polyclonal antibody targeting human SUMO2/3, validated for Western blot, IHC-paraffin, and ELISA. Detects endogenous SUMO2/3 in human, mouse, and rat samples. Suitable for studying SUMOylation and ubiquitin-like protein pathways.

[Learn More](#)

Parameter	Value
Product Name	SUMO2/3 Rabbit pAb
Remarks/Alias	HSMT3; SMT3B; SUMO3; Smt3A; SMT3H2; SUMO2/3
Species	Human
GeneID (human)	6613/6612
GeneID	6613/6612
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61956/P55854
Protein Weight	11kDa
Shipping	Ice bag