


Anti-Sumo3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55854

Item Number: CM0022206



Introduction

Rabbit polyclonal antibody against SUMO3, validated for WB, IF/ICC, and ELISA. Specifically recognizes human, mouse, and rat SUMO3, a ubiquitin-like modifier involved in nuclear transport and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	SUMO3 Rabbit pAb
Remarks/Alias	SMT3A; Smt3B; SMT3H1; SUMO-3; SUMO3
Species	Human
GeneID (Human)	6612
GeneID	6612
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human SUMO3 (NP_008867.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55854
Protein Weight	12kDa
Shipping	Ice bag