


Anti-Sumo3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P55854

Item Number: CM0026580



Introduction

Rabbit polyclonal antibody against SUMO3 (SUMO-2/SMT3A), validated for WB, IHC-P, ELISA, reacting with human, mouse, and rat. Ideal for studying ubiquitin-like modification, DNA repair, and signal transduction pathways.

[Learn More](#)

Parameter	Value
Product Name	SUMO3 Rabbit pAb
Remarks/Alias	SMT3A, Smt3B, SMT3H1, SUMO-3, SUMO3
Species	Human
GeneID (Human)	6612
GeneID	6612
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SUMO3 (NP_008867.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55854
Protein Weight	12kDa
Shipping	Ice bag