

Anti-Sumo4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Eev6

Item Number: CM0024935

Introduction

Rabbit polyclonal antibody against human SUMO4 (small ubiquitin-like modifier 4), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Suitable for studying sumoylation, oxidative stress, and immune-related pathways.

[Learn More](#)

Parameter	Value
Product Name	SUMO4 Rabbit pAb
Remarks/Alias	IDDM5; SMT3H4; SUMO-4; dj281H8.4; SUMO4
Species	Human
Gene ID (Human)	387082
Gene ID	387082
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 40-95 of human SUMO4 (Q6EEV6)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6EEV6
Protein Weight	11kDa
Shipping	Ice bag