

Anti-Sumo4 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6Eev6

Item Number: CM0008588

Introduction

Anti-SUMO4 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for SUMOylation, oxidative stress, and diabetes studies.

[Learn More](#)

Parameter	Value
Product Name	SUMO4 Rabbit mAb
Remarks/Alias	IDDM5; SMT3H4; SUMO-4; dJ281H8.4; SUMO4
Species	Human
GeneID (Human)	387082
GeneID	387082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human SUMO4 (Q6EEV6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6EEV6
Protein Weight	11kDa
Shipping	Ice bag