


Anti-Sun2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Uh99

Item Number: CM0002970



Introduction

High-quality rabbit monoclonal antibody targeting human SUN2 (UNC84B), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for LINC complex and nuclear migration studies.

[Learn More](#)

Parameter	Value
Product Name	SUN2 Rabbit mAb
Remarks/Alias	UNC84B, rab5IP, SUN2
Species	Human
Gene ID	25777 (Human)
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 618-717 of human SUN2 (Q9UH99).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UH99
Protein Weight	80kDa
Shipping	Ice bag